

# Pelican

*feminine healthcare*  
an eakin company 

Product Catalogue

Distributed by Accuratus AG

accuratus<sup>®</sup>

Distributed by Accuratus AG

**accuratus**<sup>®</sup>

Tel. +41 31 372 17 11

Fax +41 31 372 22 41

[mailbox@accuratus.ch](mailto:mailbox@accuratus.ch)

[www.accuratus.ch](http://www.accuratus.ch)

[www.subtilis.com](http://www.subtilis.com)

**Pages****Contents**

09

PELIspec® Vaginal Speculum Range

17

Scissors &amp; Forceps

22

Cervical Dilator &amp; Uterine Sound

# Welcome to Pelican Feminine Healthcare

Pelican Feminine Healthcare Ltd is a UK manufacturer of **single-use** medical instruments for **Gynaecology** and **Obstetrics**, including the market leading PELIspec<sup>®</sup> vaginal speculum.

**State-of-the-art manufacturing** facilities and active research and development highlight Pelican Feminine Healthcare's commitment to the continuous development of **innovative healthcare solutions** in Women's health.

This, combined with a positive input from specialist clinicians and other professional customers, ensures that the company remains at the forefront of providing **clinically and cost effective products.**

15,000,000  
Market-leading  
PELIspec<sup>®</sup>  
vaginal speculum  
manufactured  
since 2011

---

Pelican Feminine Healthcare  
**Production Statistics 2011 to 2014**

---

Quality  
Assurance  
Committed to  
providing quality  
products,  
manufactured to  
the highest  
standards.

Delivering market leading innovations to the highest quality standards is reinforced by **stringent monitoring** of manufacturing processes with over 15,000 quality control checks performed annually and across Pelican Feminine Healthcare's product range.

Pelican FH has invested heavily in the design and manufacture of products including the market-leading PELIspec® - **strength tested every four hours** to ensure it meets stringent product specifications. These checks are conducted at a state-of-the-art manufacturing facility in Cardiff, which since 2011 has served both UK and international markets with production of over fifteen million vaginal speculum.

The company is focused on providing quality products, manufactured to the highest standards, using **superior medical grade materials** and continuous improvement is the hallmark of Pelican FH's activities.



Pelican Feminine Healthcare are full accredited for ISO 13485:2003 and all products are CE marked in accordance with European Directive 93/42/EEC. All sterile processes are assessed for sterility by an MHRA accredited regulatory body.

Each device has a unique product identification number which can be recorded on patient notes for traceability.

All products are Latex-free, sterile, and comply with National Decontamination Guidelines.





# PELIspec® Vaginal Speculum Range



The marketing-leading PELIspec® vaginal speculum range continues to lead the way in providing a high quality single-use speculum to the healthcare market, and has been the speculum of choice for discerning medical professionals for **over 20 years**.

**Engineered and moulded in the UK** using specially chosen medical grade BPA-free (Bisphenol-A) plastic, the PELIspec® vaginal speculum range provides the ultimate combination of strength and flexibility with patient safety.

**“Fail-SAFE”** points are uniquely incorporated into the PELIspec® design, ensuring that in the event of ‘Vaginismus’ or other unanticipated stress the device will flex and break safely outside the vagina.

The PELIspec® range offers a **wide variety** of speculum solutions tailored to each procedure requirements from built-in smoke extraction and light source to the ‘Pro-Wall’ with additional wall support. The more conventional core range accommodates for every patient, available in a wide range of sizes from virgin and medium to medium-long and broad, ‘with’ and ‘without lock’.

Pelican Feminine Healthcare is committed to promoting awareness of cervical cancer and continues to donate **5p from every box of PELIspec® vaginal speculum sold**, to Jo’s Cervical Cancer Trust - the only UK charity dedicated to women and their families affected by cervical cancer and cervical abnormalities. We have now **donated over £40,000\*** to support this vital cause.


\*Figure accurate as of March 2015


Pelican are contributing 5p for every box of PELIspec sold



Registered charity in England and Wales, number: 1132642  
Registered charity in Scotland, number: SC041236

“Strong, but flexible  
design engineered  
with safety  
in mind”





# PELIspec<sup>®</sup> Vaginal Speculum

The market leading PELIspec<sup>®</sup> is a strong, sterile, disposable vaginal speculum, available in 5 sizes. Selected sizes are also available without the locking mechanism.

PELIspec<sup>®</sup> can be used for gynaecological procedures, including cervical screening, IUD fittings and vaginal examinations.

## Features

- ✓ Proven strength with built-in 'Fail-Safe' points
- ✓ Textured handles for enhanced grip
- ✓ Crystal clear for maximum visibility
- ✓ BPA and Latex-free
- ✓ Unique product identification number for traceability




Fail Safe



Manufactured in the UK



Code	Description	Units
400091	PELIspec <sup>®</sup> - Virgin with Lock	25
400105	PELIspec <sup>®</sup> - Small with Lock	25
400106	PELIspec <sup>®</sup> - Medium with Lock	25
400107	PELIspec <sup>®</sup> - Medium Long with Lock	25
400108	PELIspec <sup>®</sup> - Broad with Lock	25
400103	PELIspec <sup>®</sup> - Mixed Box with Lock	25
400104	PELIspec <sup>®</sup> - Mixed Box without Lock	25



# PELIspec<sup>®</sup> Pro-Wall

NEW

The PELIspec<sup>®</sup> Pro-Wall is the 'Professional' solution specially designed for use in cases of anterior vaginal wall laxity.

Based on the market leading PELIspec<sup>®</sup> vaginal speculum, the device features a Latex-free sheath configured to efficiently secure an unobstructed view of the cervical anatomy to ensure maximum visibility during the procedure.

## Features

- ✓ Ready applied Latex-free sheath
- ✓ Proven strength with built-in 'Fail-Safe' points
- ✓ BPA and Latex-free
- ✓ Textured handles for enhanced grip
- ✓ Unique product identification number for traceability



Manufactured  
in the UK

Code	Description	Units
400233	PELIspec <sup>®</sup> Pro-Wall Speculum - Medium	10
400235	PELIspec <sup>®</sup> Pro-Wall Speculum - Medium Long	10



# PELIspec<sup>®</sup> with Light Source

Based on the market leading PELIspec<sup>®</sup> vaginal speculum, the unique PELIspec<sup>®</sup> with Light Source features a built-in light source attachment configured to provide effective directional light for maximum visibility during gynaecological procedures.

Available in a 'Starter pack' - containing 2 x PELIspec<sup>®</sup> (1 x Medium, 1 x Medium Long), 2 x Reusable Light source leads, 2 x plastic lead covers, 1 x Re-usable Battery Pack.

## Features

- ✓ Directional light source for maximum visibility
- ✓ Compact re-usable battery pack with belt clip
- ✓ Sealed light source insertion point to ensure protection during use
- ✓ Proven strength, with 'Fail-SAFE' Points uniquely built into each PELIspec<sup>®</sup>
- ✓ Available in two sizes – Medium and Medium Long
- ✓ Reusable Leads supplied with plastic covers to reduce cross-contamination



Code	Description	Units
400207	PELIspec <sup>®</sup> with Light Source - Medium*	25
400209	PELIspec <sup>®</sup> with Light Source - Medium Long*	25
400200	PELIspec <sup>®</sup> with Light Source Starter Pack - Mixed	1
400210	Light Source Leads**	5
400220	Light Source Battery Pack	1

\*Lead covers included/ \*\*Lead covers not included



# PELIspec<sup>®</sup> with Smoke Extractor

The PELIspec<sup>®</sup> with Smoke Extractor is specially designed for use during loop diathermy procedures.

Based on the market leading PELIspec<sup>®</sup> vaginal speculum, the device features a built-in extraction tube configured to efficiently remove smoke produced during the procedure.

The PELIspec<sup>®</sup> with Smoke Extractor is deliberately manufactured in white to best reflect light and ensure maximum visibility.

## Features

- ✓ Effective smoke extraction
- ✓ Proven strength with built-in 'Fail-Safe' points
- ✓ BPA and Latex-free
- ✓ Textured handles for improved grip



Fail Safe



Manufactured in the UK

Code	Description	Units
400111	PELIspec <sup>®</sup> with Smoke Extractor - Medium	20
400112	PELIspec <sup>®</sup> with Smoke Extractor - Medium Long	20
400113	PELIspec <sup>®</sup> with Smoke Extractor - Broad	20



# PELIspec<sup>®</sup> Sims Speculum

The PELIspec<sup>®</sup> Sims Speculum is commonly used for inspection of the vaginal walls and can be used as an alternative to the traditional Cusco design.

The design and shape of the PELIspec<sup>®</sup> Sims Speculum is characterised by its two right-angled bends and differing sized ends.

## Features

- ✓ Crystal clear ensuring greater visibility of the vaginal wall
- ✓ Smooth rounded ends for patient comfort
- ✓ Mirrors design of metal equivalent
- ✓ Easy central grip handle
- ✓ Two differing end sizes for convenience



Manufactured  
in the UK

Code	Description	Units
400131	PELIspec <sup>®</sup> Sims Speculum - Medium	25
400132	PELIspec <sup>®</sup> Sims Speculum - Medium	3
400135	PELIspec <sup>®</sup> Sims Speculum - Large	10



# Scissors & Forceps





# Pelican Spongeholder

The Pelican Spongeholder is designed to ensure visibility is not impeded when used in conjunction with a PELISpec® vaginal speculum.

The robust design gives the user full control of the device and the grooved tips ensure a secure grip of swabs.

## Features

- ✓ Grooved tip for maximum grip
- ✓ Enhanced secure locking mechanism
- ✓ Specifically designed for use with swabs
- ✓ Cost effective alternative to metal with parallel performance



Manufactured  
in the UK

Code	Description	Units
400300	Pelican Spongeholder	50
400302	Pelican Spongeholder	5



# Pelican Vulsellum

The Pelican Vulsellum is characterised by a 3:4 toothed grip specially designed to offer a gentle but secure grip on the cervix.

The Pelican Vulsellum has an angled head to maximise access to the cervix during the procedure and improve visibility.

Provides an alternative choice for practitioners who prefer a more evenly spread grip over a singled toothed (Tenaculum) grip.

## Features

- ✓ Enhanced secure locking mechanism
- ✓ Specifically designed tooth combination to minimise discomfort and reduce trauma
- ✓ Angled tips allow greater visibility
- ✓ Cost effective alternative to metal with parallel performance



Manufactured  
in the UK

Code	Description	Units
400335	Pelican Vulsellum	50
400336	Pelican Vulsellum	5



# Pelican Tenaculum

The Pelican Tenaculum features a precise single toothed grip with secure locking mechanism, ideal for holding the anterior lip of the cervix during procedures.

The Tenaculum provides an alternative choice to practitioners who prefer a single toothed grip over a more evenly spread grip.

## Features

- ✓ Precise single toothed grip
- ✓ Secure locking mechanism
- ✓ Cost effective alternative to metal with parallel performance
- ✓ Specifically designed to allow greater visibility of the cervix



Manufactured  
in the UK

Code	Description	Units
400330	Pelican Tenaculum	50
400331	Pelican Tenaculum	5

# Pelican Coil Remover

The Pelican Coil Remover is specifically designed to remove intra-uterine devices and is characterised by an angled tip with a 2:1 tooth combination.

The robust design gives the user full control of the instrument and the “V” designed tip ensures a secure grasp of the coil thread to facilitate retrieval of the intra-uterine device. The tip size has been reduced and angled for optimal visibility during the procedure.

## Features

- ✓ Secure locking mechanism
- ✓ Angled tip to optimise visibility
- ✓ Specially designed tooth combination facilitates easy IUD retrieval



Code	Description	Units
400317	Pelican Coil Remover	50
400318	Pelican Coil Remover	5

# Cervical Dilator & Uterine Sound





# Pelican Uterine Sound-Dilator

NEW

The Uterine Sound-Dilator (“USD”) is a strong but flexible, single-use, double-ended combined device designed for gentle dilation of the cervix and to assess the depth of the uterus.

The smooth and rounded ends are specially designed to allow gradual dilation of the cervix, whilst reducing risk of perforation.

This ‘2-in-1’ instrument, provides added convenience when fitting intra-uterine devices

## Features


- ✓ Strong but flexible moulded plastic
- ✓ Ideal for coil-fitting procedures
- ✓ Smooth, rounded tips for reduced risk of perforation
- ✓ Clear and precise 1cm depth indicators
- ✓ 27.7cm long with a maximum diameter of 5.5mm



2 in 1



Manufactured  
in the UK



# Pelican Dilator

The Pelican Dilator is designed to gently dilate the cervical 'os' to facilitate access to the uterus, and is ideal for preparing the cervix for the introduction of an intra-uterine device.

The Pelican Dilator is tapered to allow gradual dilation in 1mm increments from 2mm to 6mm (using two devices).

## Features

- ✓ Strong but flexible moulded plastic
- ✓ Central thumb hold for ease of use
- ✓ Tapered ends allow gradual cervical dilation
- ✓ Smooth, rounded tips for reduced risk of perforation



Manufactured  
in the UK

Code	Description	Units
400342	Pelican Dilator 3/4mm	100
400343	Pelican Dilator 3/4mm	10
400344	Pelican Dilator 5/6mm	100
400345	Pelican Dilator 5/6mm	10



# Pelican Uterine Sound

The Pelican Uterine Sound is specifically designed to confirm uterine depth prior to the insertion of an intra-uterine device. The rounded end ensures a gentle introduction of the device through the cervical 'os' into the uterine cavity.

The Pelican Uterine Sound is a slender, flexible gauge with 1cm indicators to allow immediate visual confirmation of uterine depth.

## Features

- ✓ Notches allow measurement of the uterine depth
- ✓ White in colour to aid visualisation
- ✓ Mirrors the metal equivalent for ease of use
- ✓ Increased patient comfort due to design
- ✓ Unique product identification number for traceability



Code	Description	Units
400310	Pelican Uterine Sound	100



# PELIhook Amniotic Membrane Perforator “Amnihook”

The PELIhook Amniotic Membrane Perforator is designed to safely rupture the amniotic membrane whilst minimising the risk of injury to the foetus.

The sharp hook cleanly snags the membrane, and the rounded tip ensures that the hook does not inadvertently come into contact with maternal tissue.

## Features

- ✓ The curved PELIhook is anatomically shaped to allow ease of use
- ✓ Designed to safely rupture the membrane
- ✓ Textured thumb grip on the handle
- ✓ Option of a straight PELIhook for optimal control and comfort
- ✓ Sterile and single-use
- ✓ Latex-free



Code	Description	Units
SH2	PELIhook Amniotic Membrane Perforator - Straight	100
KH2	PELIhook Amniotic Membrane Perforator - Curved	100



# Pelican FETOfit®

NEW

FETOfit® is a disposable device specifically designed to secure both 'TOCO' and 'CTG' transducers quickly and easily for strapless foetal monitoring.

Unlike standard circumferential straps, FETOfit's hypoallergenic adhesive pads adhere directly to the front of the abdomen.

The flexible material moulds to body contours, offering a secure fit whilst promoting maternal comfort and mobility whilst facilitating access to the back in the event of an epidural.

## Features

- ✓ Universal attachment to suit all transducers
- ✓ Quick and easy to re-position
- ✓ Easy access for epidural
- ✓ Disposable device for single-person use
- ✓ Secure fit
- ✓ Safeguards patient comfort & mobility



Manufactured  
in the UK

Code	Description	Units
400625	FETOfit®	25



# PELIhook Amniotic Membrane Perforator “Amnihook”

The PELIhook Amniotic Membrane Perforator is designed to safely rupture the amniotic membrane whilst minimising the risk of injury to the foetus.

The sharp hook cleanly snags the membrane, and the rounded tip ensures that the hook does not inadvertently come into contact with maternal tissue.

## Features

- ✓ The curved PELIhook is anatomically shaped to allow ease of use
- ✓ Designed to safely rupture the membrane
- ✓ Textured thumb grip on the handle
- ✓ Option of a straight PELIhook for optimal control and comfort
- ✓ Sterile and single-use
- ✓ Latex-free



Code	Description	Units
SH2	PELIhook Amniotic Membrane Perforator - Straight	100
KH2	PELIhook Amniotic Membrane Perforator - Curved	100

Distributed by Accuratus AG

**accuratus**<sup>®</sup>

Tel. +41 31 372 17 11

Fax +41 31 372 22 41

[mailbox@accuratus.ch](mailto:mailbox@accuratus.ch)

[www.accuratus.ch](http://www.accuratus.ch)

[www.subtilis.com](http://www.subtilis.com)

**Pelican**  
*feminine healthcare*  
an eakin company 

Manufacturer:

**Pelican Feminine Healthcare Ltd.**  
Greypoint, Cardiff Business Park, Cardiff CF14 5WF

Distributed by Accuratus AG

**accuratus**<sup>®</sup>

Tel. +41 31 372 17 11

Fax +41 31 372 22 41

[mailbox@accuratus.ch](mailto:mailbox@accuratus.ch)

[www.accuratus.ch](http://www.accuratus.ch)

[www.subtilis.com](http://www.subtilis.com)